

<b>Model</b>	<b>685BF1</b>	<b>785BF1</b>	<b>695BF1</b>	<b>795BF1</b>
<b>Description</b>	8LPM 115V	8LPM 230V	10LPM 115V	10LPM 230V
<b>Frequency</b>	60Hz	50/60Hz	60Hz	50Hz
<b>Average Power</b>	500 Watts	490 Watts	700 Watts	600 Watts
<b>Protection Class</b>	Class II			
<b>Mains Protection</b>	10A or 13A	5A or 7A	10A or 13A	5A or 7A
<b>Liter Flow</b>	2 to 8 LPM		2 to 10 LPM	
<b>Outlet Pressure</b>	15 PSIG		20 PSIG	
<b>Average Oxygen Content</b>	At 2 LPM > 90%		At 2 LPM > 90%	
<b>Average Oxygen Content</b>	At 8 LPM 87% to 95.5%		At 10 LPM 87% to 95.5%	
<b>Low O<sub>2</sub> Liter Flow</b>	80 to 90 LPM		90 to 100 LPM	
<b>Average Low O<sub>2</sub> Content</b>	At 2 LPM on Unit 15% to 18%		At 2 LPM on Unit 15% to 18%	
<b>Average Low O<sub>2</sub> Content</b>	At 8 LPM on Unit 10% to 14%		At 10 LPM on Unit 10% to 14%	
<b>Dimensions (L x W x H)</b>	39 x 40 x 71 cm (15 x 15.5 x 28 in.)			
<b>Weight</b>	24.5-26 kg (54-58 lbs.)			